



RUBELLA IgG ANTIBODIES

Date: February 12, 2008

Subject: Method Change for Rubella IgG Testing

To better serve you, the EMH Reference Laboratory has moved Rubella IgG testing to an automated method. This change will improve turn around time for this test. The specimen requirements will not change, however please note the change in reference range in the table below. Please update your Client Users Manual with these changes.

Test Name	Method	Specimen Requirements	Reference Range IU/mL	TAT	Order code	CPT
Rubella IgG	Chemiluminescent Immunoassay	1 SST 0.5 ml serum Refrigerate (>48 hrs, freeze)	<10 Non-immune 10-14.9 Equivocal >15 Immune	24 hrs daily	RUBG2	86762

For more information regarding this notice, please call Client Services at 866-941-4542.

jmh Feb-08